




Can you align a tibia using  
finesse, and not force?



 **smith&nephew**  
**TRIGEN<sup>◇</sup>**  
**META-NAIL<sup>◇</sup>**  
Semi-extended Instrument Set



Through constant innovation, Smith & Nephew is the **first company to fully support** an alternative technique to the standard nailing of tibia fractures.



Easier alignment

Simplified procedure


Enhanced fixation

The new TRIGEN<sup>®</sup> Semi-extended Instrumentation facilitates intramedullary nailing procedures of the tibia with **less flexion of the knee joint**. Reduced flexion leads to less pull from the quadriceps mechanism, helping to **avoid fracture malreductions and misalignment**.

The semi-extended technique offers easier reductions, easier imaging and less angulation of the proximal segment while offering the **enhanced, multiplanar fixation** synonymous with the META-NAIL<sup>®</sup> Tibial Nailing System.




## Challenges with tibial nailing



133 tibia nailing cases showed  
58% of proximal tibia fractures  
were malaligned (>5° angulation)

# Easier alignment



Stresses from angle of flexion in traditional technique

When dealing with proximal tibia fractures, a common issue is forced angulation. The semi-extended position allows for **neutralization of the pulling forces of the quadriceps tendon.**<sup>2</sup>

Without these pulling forces the fracture **naturally reduces**, leading to **less malalignment.**



Reducer used in extended position allows for easier control of the distal fragment.


# Simplified procedure

**C-Arm use becomes easier** for M/L & A/P views due to leg position.

Leg positioning does not need to change for freehand distal locking.

**Reaming and guide wire placement are easier** since they do not have to be raised in the air over the hyper-flexed knee.

There is no need for triangles or other aids to force the position of the leg. In the semi-extended position **only bumps are needed for the procedure.**




# Enhanced fixation – Proven META-NAIL<sup>◇</sup> technology

The Semi-extended Instrumentation still utilizes all the benefits of the TRIGEN<sup>◇</sup> META-NAIL System.

Tissue protecting instruments built specifically for this approach allow surgeons to operate in the joint space with confidence and accuracy.

The **multiplanar screw configuration** of the META-NAIL Tibial Nail continues to offer a **stable, locked construct**.




With the semi-extended instruments, the surgeon is still able to use the **blocking screw attachments** to create an artificial diaphyseal canal and **help prevent nail translation** as the patient begins to move and weight bear.




Adjustable cartridge allows canal to be lengthened

The threaded holes and polyethylene bushing **enhance screw fixation** and help **prevent screw back out and nail translation**.



As the proximal screw is moved down the dynamization slot through the guide, the distal segment of the tibia is pulled up to reduce the fracture.



The META-NAIL system allows for up to 7mm of compression while keeping the guide attached. This leads to **more control of the reduction and easier, reproducible outcomes**.

The META-NAIL Tibial Nail continues to offer the highest proximal bend on the market, allowing for easier insertion into the canal space with less risk of dislocation.

think.  
again

## References

1. Freedman EL, Johnson EE. Radiographic analysis of tibial fractures malalignment following intramedullary nailing. *Clin Orthop Relat Res.* 1995:25-33.
2. Tornetta P, 3rd, Collins E. Semi-extended position of intramedullary nailing of the proximal tibia. *Clin Orthop Relat Res.* 1996:185-189.

### **Orthopaedic Reconstruction & Trauma**

Smith & Nephew, Inc.  
1450 Brooks Road  
Memphis, TN 38116  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Telephone: 1-901-396-2121  
Information: 1-800-821-5700  
Orders/Inquiries: 1-800-238-7538