



 **smith&nephew**
TAYLOR SPATIAL
FRAME[◊]
External Fixator

The revolutionary yet simple solution
for trauma and deformity correction.

The TAYLOR SPATIAL FRAME Fixator is the world's most advanced external
fixation system, combining ease of application, proven stability and
effectiveness with the unmatched precision of computer-calculated corrections.



Versatile

Its unique six-strut configuration allows one frame to correct deformities, and the new Fast Fx struts allow the frame to be used for acute trauma applications. The TAYLOR SPATIAL FRAME Fixator can correct multiplanar angular and transitional deformities simultaneously.

Easy to use

The TAYLOR SPATIAL FRAME Fixator is a true computer-assisted trauma product, driven by the world's most powerful operative planning software. The software is simple and easy to use, and a Help Desk is available 24 hours a day, seven days a week. Once data is entered, calculations are completed in seconds, providing the exact specifications for strut adjustments to obtain precise alignment.

Precise fracture alignment

The TAYLOR SPATIAL FRAME System gives you power of precision in treating deformities and acute trauma. The unique frame construct plus the software's Total Residual method allows you to achieve the level of precision alignment you desire. Your lateral data input is summarized on a screen and another view displays the frame and bone status so you can easily and quickly check your patient's progress. The program monitors the fixator and identifies needed adjustments to the struts. When acute fracture reductions cannot be accomplished, the TAYLOR SPATIAL FRAME Fixator allows you to make the reduction in stages, so that precise fracture alignment is eventually achieved.

The TAYLOR SPATIAL FRAME Fixator. The world's most advanced external fixator.

Log on to www.spatialframe.com to register for an account. Call Smith & Nephew Medical Education at 1-800-544-9672 for information on the TAYLOR SPATIAL FRAME Fixator product training courses.

TAYLOR SPATIAL FRAME™ smith&nephew

Home Cases Utilities Literature

NEW! The Taylor Spatial Frame Fixator VERSION 2.1™

The simple solution for treating acute fractures and complex deformities is finally within your reach. The Taylor Spatial Frame Fixator.

At the heart of the Taylor Spatial Frame System is an algorithm created to perform the calculations needed to accurately manipulate the Taylor Spatial Frame, a 6 degrees-of-freedom external fixator. This algorithm, providing you with visual feedback, is available to you via this website to enable an efficient and successful correction... Getting it right the first time! This easy to use Web-based Application guides you through every step of the way.

SPATIAL NEWS

Version 2.1 is now available! Click on "What's New" to see the changes.

Please join us at our next Spatial Frame course February 20-22 in Snowmass, Colorado. For course information, please call Bonnie Muse at 901-399-5202 or e-mail at bonnie.muse@smithnephew.com

Secure Customer Login

User Name

Password

[Request an Account](#) [Forgot User Name or Password? Click here](#)

TSF Web Application Version 2.1.0
Taylor Spatial Frame™ property of [Smith & Nephew](#). All Rights Reserved.
Application Designed and Managed by [NetIDEAS, Inc.](#)

TAYLOR SPATIAL FRAME™ smith&nephew

Home Cases Utilities Literature User Profile Logout Help

File Case Info Define Deformity Select Frame Mount Frame Initial Frame Final Frame Structure at Risk Prescription Report

Initial Settings for Total Residual Operative Mode

Case Name: GEA HTO Surgery 8/12/02

Strut 1 (mm) (Red)	Strut 2 (mm) (Orange)	Strut 3 (mm) (Yellow)	Strut 4 (mm) (Green)	Strut 5 (mm) (Blue)	Strut 6 (mm) (Violet)
150	145	140	135	160	160

Right AP View Right Lateral View Right Axial View

Lateral Medial Posterior Anterior Medial Lateral

Deformity Parameters

AP View Angulation: 23.0° Varus
Lateral View Angulation: 10.0° Apex Anterior
Axial View Angulation: 15.0° Internal

AP View Translation: 10.0 mm Medial
Lateral View Translation: 5.0 mm Posterior
Axial Translation: 5.0 mm Short

Mounting Parameters

AP View Frame Offset: 0.0 mm
Lateral View Frame Offset: 20.0 mm Posterior to Origin
Rotary Frame Angle: 0.0°
Axial Frame Offset: 30.0 mm Proximal to Origin

Clicking on graphic will enlarge ?

TSF Web Application Version 2.1.0
Taylor Spatial Frame™ property of [Smith & Nephew](#). All Rights Reserved.
Application Designed and Managed by [NetIDEAS, Inc.](#)

TAYLOR SPATIAL FRAME™ smith&nephew

Home Cases Utilities Literature User Profile Logout Help

File Case Info Define Deformity Select Frame Mount Frame Initial Frame Final Frame Structure at Risk Prescription Report

Final Settings for Total Residual Operative Mode

Case Name: GEA HTO Surgery 8/12/02

Strut 1 (mm) (Red)	Strut 2 (mm) (Orange)	Strut 3 (mm) (Yellow)	Strut 4 (mm) (Green)	Strut 5 (mm) (Blue)	Strut 6 (mm) (Violet)
157	131	198	166	152	121

Right AP View Right Lateral View Right Axial View

Lateral Medial Posterior Anterior Medial Lateral

Final Deformity Parameters

AP View Angulation: 0.0°
Lateral View Angulation: 0.0°
Axial View Angulation: 0.0°

AP View Translation: 0.0 mm
Lateral View Translation: 0.0 mm
Axial Translation: 0.0 mm

Mounting Parameters

AP View Frame Offset: 0.0 mm
Lateral View Frame Offset: 20.0 mm Posterior to Origin
Rotary Frame Angle: 0.0°
Axial Frame Offset: 30.0 mm Proximal to Origin

Clicking on graphic will enlarge ?

TSF Web Application Version 2.1.0
Taylor Spatial Frame™ property of [Smith & Nephew](#). All Rights Reserved.
Application Designed and Managed by [NetIDEAS, Inc.](#)