


# Need a quick fix?

Introducing the TAYLOR SPATIAL FRAME  
FAST FX<sup>◊</sup> Strut.

 **smith&nephew**  
**TAYLOR SPATIAL**  
**FRAME<sup>◊</sup>**  
External Fixator



# A quick fix. The TAYLOR SPATIAL FRAME FAST FX Strut.

The TAYLOR SPATIAL FRAME FAST FX struts make fracture reduction faster and easier than ever before. In the acute reduction mode, the struts allow the frame to be used as a fracture reduction device. Attach the frame to the bone. Unlock the struts. Reduce the fracture. Lock the struts. You're done. It's just that simple.

Still need more correction? Just run a residual. In the gradual mode, the FAST FX struts allow for fine tuning until the correction is complete. FAST FX delivers the perfect corrections you expect from the TAYLOR SPATIAL FRAME.

## Features

Locking adjustment knob can be disengaged for acute adjustments

Gradually adjustable

Jam nut seats against the adjustment knob

ID band clips on adjustment knob

Available in four sizes

Large print on strut body and ID bands

## Benefits

Faster and easier adjustments in all six planes; when all six struts are disengaged, the frame can be used as a fracture reduction device

Adjustments can be fine-tuned as precisely as you want

Maintains strut adjustments by eliminating strut migration

Prevents accidental disengagement of adjustment knob

Accommodates a broad range of anatomy and strut adjustments

## Catalog Information

### FAST FX Strut

Cat. No.	Description	Adjustment Range
7107-0705	X-Short	91-121mm
7107-0710	Short	116-152mm
7107-0720	Medium	143-205mm
7107-0730	Long	195-311mm



### Standard Strut

Cat. No.	Description	Adjustment Range
7107-0205	X-Short	75-96mm
7107-0210	Short	90-125mm
7107-0220	Medium	116-178mm
7107-0230	Long	169-283mm



### ID Bands

Cat. No.	Description
7107-0340	ID Band Kit



### ID Bands

Cat. No.	Description
7107-0320	ID Band Kit